

SELECTION OF 2 SOLVENTS



RINING AND CONDITIONING OF HPLC COLUMNS

INVER-CIL®

APPLICATION

Immediate automatic change from solvent A to solvent B

- by push button
- by electric signal given by a pump or integrator



ADVANTAGES

- Works in your absence.
- Avoids repetitive manipulations of harmful products
- Low cost
- Small dimensions
- No compressed air required

SPECIFICATIONS INVER-CIL®	
Power supply	230V
Communications	Manual, Remote
Dimensions	14 x 7 x 4 cm
Weight	400 g

MODEL	REFERENCE
INVER-CIL®, 230V 2 solvent selector	99S1

amchro GmbH

Chromatographie- und Laborzubehör

Hugo-von-Eltz-Str. 22
D-65795 Hattersheim

TEL : (0800) 26 24 768 / FAX : (0800) 26 24 763
E-mail : info@amchro.com